HELIO
Olympia[®]
Kinesiology Tape

HELIO
Olympia
Kinesiology Tape

Support joint and muscle.

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Best for injured **muscles** or
joints, swelling and pain relief.

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Expanding the space between **skin**
and **muscles** to release **pressure**
and **improve** fluid circulation.

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Relief **tender tension**.
Improved poor **posture** and
breathing.

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光沢のある新素材
Fabric reflecting
under the lights.

Point2


日本製粘着のり Japanese made top grade acrylic adhesive.

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皮膚からの水分や汗が通り抜ける
Superb water and sweat resistant characteristics.

〈粘着面〉 Underneath

粘着面が皮膚を引き上げる
Expanding the space between skin and muscles.

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How to Use the PRE-CUT type (25cm/cut)

PRE-CUT type contains a cut in 25cm between, it can be used immediately without scissors.

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4 Gently pull both sides to separate the tape.
軽く引っ張るとテープが切れます。

How to Cut the Tape in Half

Sometimes depend on the taping position you may need to cut the tape in half.

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Without Tension to the tape at the Beginning and the End when taping

Without pulling tension in the beginning and the end the tape is less likely to be wiping off.

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04 Reduce Fatigue at Front of the Shin Bone
スネの前側の疲労軽減
→ P.39

ANKLE&FOOT — 足首・足 —

01 Improved Stability of Ankle
足首の安定性を向上
→ P.40

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
SHOULDER 01 —肩関節—

Improve Shoulder Movement

Gently support shoulder area improved stability

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1
Begin the tape at the 1/2 of the upper arm tape along from the back of the shoulder following the muscle line.

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2
Start 2nd tape as the same location as the 1st, covered the edge of 1st tape following the front muscle line with light tension ended the tape on the shoulder line.

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Complete image

SHOULDER 02 —肩関節—

Reducing Pain from Back of the Shoulder

Reducing pain or discomfort from back of the shoulder.

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
1
Begin by placing the tape at the front of the shoulder, applying the tape cover trapezius central area with light tension.

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2
Use half of the Per-Cut tape applying the tape in the back of the shoulder with light tension from down to up.

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3
Use the other half of the tape repeat step 2 parallel by cover 1/5 of step 2.

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Complete image

SHOULDER 03 —肩関節—

Increase Shoulder Joint Stability

Increase stability and support of shoulder joint.

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1
Use 2 pieces half of the Per-cut applying the tape cross from down to up with light tension.

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2
Repeat step 1, applying from different direction from down to up with light tension.

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
ELBOW 01 —肘—

Increase Elbow Stability

Reduce pain while extend elbow, improved elbow stability.

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
1
Position the elbow 90 degrees start the tape from the middle of forearm around the outer side of the elbow ended the tape at the upper arm.

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2
Keep the same position repeat step 1, applying the tape at the inner of the arm.

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ELBOW 02 —肘—

Inner of the Elbow Support

Reduce pain from inner elbow for overuse muscle and tender from sports.

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Complete image

1

Use half of the Pre-cut by position elbow 90 degrees. Apply tape on the elbow cover the ligament to the inner of the upper arm.

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2

Use the other half of the Per-cut tape keep the same position applying the tape from inner of the elbow cross over the ligament from up to done.

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3

Start the tape at the outer of the elbow crossing the ligament ended the tape at inner of the forearm close to the wrist.

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
ELBOW 03 —肘—

Outer Elbow Support (Tennis Elbow)

Reduce pain from outer elbow for overuse muscle or tender from sports.

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Complete image

1

Use half of the Per-cut Position elbow 90 degrees. Apply the tape on the outer of the elbow.

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2

Use the other half of the Per-cut tape from up to down cross over the ligament.

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3

Start the tape at the outer of the elbow crossing the ligament ended the tape at the outer of the forearm close to the wrist.

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
ELBOW 04 —肘—

Biceps Support

Supporting biceps muscle when need heavy lifting.

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Complete image

1

Start the tape on the inner elbow hinge joint area applying the tape towards to the shoulder.

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2

Ended tape at the shoulder with mid tension.

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ELBOW 05 —肘—

Triceps Support

Supporting triceps and to balance with biceps while needed.

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Complete image

1

Begin by placing the tape on the outer elbow area.

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2

Applying the tape following by the triceps muscle ended at the lower shoulder area with mid tension.

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WRIST&FINGER 01 —手首・指—

Wrist Stability

Increase performance and stability of the wrist.

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Complete image

1
Use half of the Per-cut by placing the tape at the inner of wrist.

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2
Use the other half of the Per-cut by placing the tape on the outer of the wrist.

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3
Start the tape on the hand applying the tape towards to the elbow.

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
Complete image

1
Cut Per-Cut tape parallel half start the tape from edge of the nail following taping towards to the elbow.

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2
Start the 2nd half tape at the first thumb joint applying the tape straight down towards to the elbow.

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3
Use half of the Per-cut by placing the tape on the outer of the wrist.

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WRIST&FINGER 03 —手首・指—

Jammed Index Finger

Taping for first joint of jammed finger to stabilized and pain reduce.

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Complete image

1

Cut Per-Cut tape parallel half start the tape from the edge of the nail applying the tape straight down towards to the elbow.

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2

Cut tape in half from step 1 applying the tape on the first joint of index finger.

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
Lower Back Support

For daily work out or heavy lifting to support lower back area.

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Complete image

1

Start the tape transverse direction on the illocostalls lumborum area.

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2

Repeat step 1 by covering proximally 1cm.

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3

Cut Per-cut in half placing tapes from both side of the waist towards to the center of the back.

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Pinpoint Pain

For area pinpoint pain directly application to protect injure area for early recovery.

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1
Make sure the pain point. Cut Per-cut in half applying tape cross over.

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2
Repeat step 1 from different direction.

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Support Back Movement or Reduce Pain

For stiff muscle from lack of exercise or from stress.

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1
Start the tape following by the side of the spine make sure to inhale before applying the tape.

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2
Repeat step 1 on the other side of the spine.

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Improve Stability of the Twisting Motion

Support of twisting motion, improve stability of the lower abdomen.

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Ready for tape < 準備するテープ >

1 Start the tape on the side of the body (rib cage) twist body towards to the back while applying ended the tape by the side of the belly button.

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2 Start at the side of the step one start\ t point repeat as step 1 ended the tape on the side of the belly button.

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
Discomfort from Outer Hip

Relief outer hip discomfort or pain from walking, jogging.

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
Ready for tape < 準備するテープ >

1 Cut Per-cut in half start the tape at upper thigh from down to up, cross over the 2nd tape.

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2 Repeat step 1 from different direction.

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Hip Support

Relief muscle pain from hip or hip joint.

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
1
Start tape at the pelvis applying the tape diagonal towards to the outer thigh.

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Quadriceps Muscles Support

Could be use while doing sports, up and down stairs or hills.

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1
Begin at setting position. Start tape at the upper knee inner of the quadriceps.

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2
Repeat step on the other side of the outer quadriceps.

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Adductor Muscle Support

Discomfort during operations such as raid dash and stop.

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Ready for tape < 準備するテープ >

Complete image

1
Begin at setting position.
Start tape at the side of knee
applying the tape from inner of
the thigh.

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2
Repeat as step 1 above step 1
cover about 1cm.

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Hamstring Muscles Support

To support and stabilized hamstring muscles.

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Ready for tape < 準備するテープ >

Complete image

1
Have the lunge position ready.
Start the tape behind the knee
towards to the tibia.

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2
Keep the same position apply
the tape on the other side
towards to the tibia.

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KNEE 01 —膝—

Support of the Knee

To support the muscle and joint of the knee and improved the stability of the patella.

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Complete image

1
Begin placing the tape from bottom of the patella with light tension following the vastus lateralis towards the thigh with no tension.

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KNEE 02 —膝—

Reduce Pain of Jumper's Knee

Reduce pain and support while doing intense training or sports.

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
Complete image

1
Cut Per-Cut in half placing the tape right under the knee.

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2
Start the tape at the middle under the patella following the vastus lateralis towards the thigh.

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KNEE 03 — 膝 —

Support Lateralis of Knee

Support and stabilized overuse lateralis of knee.

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
1
Cut Per-cut in half. Position the knee in 45 degrees start the tape on the side of the patella from down to up cover the tendon ligament.

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2
Use the other half of the Per-Cut tape cross over the first tape from down to up to cover the tendon ligament.

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3
Start the tape in between Step 1 and 2, tape over the tender ligament following the lateralis thigh.

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KNEE 04 — 膝 —

Support Medialis of Knee

Support and reduce pain of medialis of knee.

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

1
Cut Per-cut in half. Position the knee in 45 degrees start the tape on the side of the patella cover the tendon ligament.

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3
Start the tape in between step 1 and 2, tape over the tendon ligament following the medialis thigh.

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Complete image

KNEE 05 — 膝 —

Pain Reduce for OSD (Osgood-Schlatter Disease)

Reduce pain from patella ligament below the knee.

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1
Cut Per-Cut in half. Start the tape diagonally on the top of the tibial tuberosity.

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2
Use the other half of the Per-Cut tape, tape cross over the step 1.

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Complete image


KNEE 06 — 膝 —

Support for Bending and Stretching of the Knee

Reduce pain and discomfort from knee and hamstring muscles

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1
Start tape from lateralis side of the knee towards to the hip.

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2
Repeat step 1 ended the tape towards to the tibia.

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Relive Calf Fatigue

Reduce the fatigue and stiffness of the calf muscles.

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1
Start the tape at the outer of the Achilles tendon following the calf muscle line.

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2
Repeat the step on the inner side of the calf muscle.

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Reduce pain from Achilles Tendon

Reduce discomfort and pain of the Achilles tendon due to overuse and fatigue.

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1
Position with heel up apply the tape on the heel cover the Achilles tendon with no tension.

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2
Cut Per-Cut in half placing the tape on the Calcaneal tuberosity.

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3
Use the other half of the Per-cut placing above step on the Achilles tendon.

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Shin Pain Relief


Reduce pain from shin bone mostly from growth phase.

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Ready for tape < 準備するテープ >

Complete image

1
Start the tape at the inner of the ankle tape along the inside of the shin bone.

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2
Cut Per-Cut in half placing the tape directly on the pain point.

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3
Use the other half the tape placing above the step 2 by covering approximately 2cm.

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Reduce Fatigue at Front of the Shin Bone

Reduce pain or discomfort on the front side of the shin bone due to ankle fatigue.

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Ready for tape < 準備するテープ >

Complete image

1
Start the tape at the front of the ankle following the shin bone.

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2
Cut Per-Cut in half placing the tape cross step 1 in the middle of the shin bone.

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3
Use the other half of the Per-Cut place the tape above the ankle cross the step 1.

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ANKLE&FOOT 01 — 足首・足 —

Improved Stability of Ankle

Improve stability of the ankle against such as foot pain or from sprain.

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
Complete image

1
Start placing the center the tape on top of the heel, ended both end at the inner and outer of the ankle.

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2
Start the tape at the bottom of the thumb passing through achillies ended at the other side of the foot.

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3
Placing the tape at the bottom of the foot closed to the toes by passing over the heel toward to the Achilles tendon without tension.

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
ANKLE&FOOT 02 — 足首・足 —

Reduce Inner Ankle Pain

Reduce pain from inner of the ankle due to from fatigue or sprain.

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Complete image

1
Start the tape at the bottom of the foot covering the arch toward to the Achilles Tendon.

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2
Repeat as step 1 placing the tape below by covering approximately 2cm.

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ANKLE&FOOT 03 — 足首・足 —

Pain from the Thumb

Pain from thumb due to overuse or fatigue

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1

Cut Per-Cut in half and cut another parallel half and placing the tape around thumb.

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2

Cut Pre-Cut in Parallel half start the tape at the bottom end of the thumb toward to the Achilles tendon.

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ANKLE&FOOT 04 — 足首・足 —

Reduce Pain from Hallux Valgus

Reduce pain and discomfort from the bottom of the thumb

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1

Cut Per-Cut in parallel half, starting the tape at the inner of the thumb toward to the back of the heel.

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2

Repeat as step 1 by placing the tape above.

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3

Cut Per-Cut tape in half, placing both tape at the top of the foot by covering the bottom of the thumb joint.

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Complete image


ANKLE&FOOT 05 — 足首・足 —

Reduce Fatigue from Bottom of the Foot

Reduce pain or discomfort from such as Plantar tendon meningitides.

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2
Start the tape at the inner side of the foot by passing the calcaneal tuberosity and the arch ended the tape at the same side as started.

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3
Start the tape at the outer side of the foot passing by the calcaneal tuberosity.

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ANKLE&FOOT 06 — 足首・足 —

Reduce Pain from Bottom of the Toes

Reduce discomfort from bottom of the toes while walking.

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Complete image

1
Cut Per-cut in Half Placing the tape at the inner side of the foot by crossing the bottom of the foot ended at the other side.

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2
Repeat Step 1 by placing the below it and covering approximately half of the tape 1.

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3
Start the 3rd tape at the middle of the tape 1 and 2 toward to the Achilles tendon without tension.

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Rounded corner 25cm cutting Pre-Cut type, Easy tear.
Continue Roll type

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5m roll type with 5cm width / 17 pieces pre-cut type / 31m
bulk roll type also has 8 colors available / 4 difference width
to choose / Compact type to commercial bulk roll type.

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4 different widths are available
for bulk pack in beige color.

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816 Charcot Ave, San Jose, CA95131, U.S.A
ITO CO., LTD.