

GENERAL CATALOGUE

for

# ACUPUNCTURE

**ITO** PHYSIOTHERAPY  
& REHABILITATION

# INDEX

MOXA ----- 1

NON MAGNETIC PRESS PELLETS ----- 9

MAGNETIC PRESS PELLETS----- 11

ACUPUNCTURE DEVICES----- 12

KINESIO ELASTIC TAPE----- 14



# MOXA

## MOXA BOX



### MOGUSAN'S MOXA BOX

No smoke & No odor Dedicated carbonized moxa includes 50 pieces.

- Achieved a surprising ease of ignition
- Set item Mogusan's Moxa Box Dedicated carbonized moxa

Using only moxa as the raw material, it ensures an ignition quality reminiscent of pure moxa. Since there is no smoke during combustion, you can use it anywhere without restriction.




Product movie



### Characteristics:

- Rotating vent hole adjustment plate for temperature control
- A mesh basket for easy ash disposal
- Connectivity feature of 2 boxes for broader application



Type	pcs/box	Order Code
Set item: Mogusan's Moxa Box	50	8009377
Refill item: Dedicated carbonized moxa	100	8009378

### Specifications

Dimensions: 110 (W) × 110 (D) × 100 (H) mm  
Weight: 300 g  
Duration of warmth: Approx. 20 minutes  
(for 1 pc of dedicated carbonized moxa)  
Material: [Main unit] Hinoki wood and ceramic  
[Mesh basket] stainless steel

# MOXA

## SMOKELESS, FRUIT FRAGRANCE STICK-ON MOXA



pcs/box	Order Code
300	8009368
1,000	8009369




### CHOSEI-KYU MOXA SMOKELESS & FRUIT FRAGRANCE

Generate slight fruit fragrance and little smoke.  
During use, no ventilation is required, and curtains and clothes are free from characteristic moxa odor.

#### Characteristics:

- Little smoke: Smoke will disappear after it is lit
- Fruit fragrance: No characteristic moxa odor
- Stable heat temperature: Besides moxa, charcoal is used as heat source.

## STICK-ON MOXA



### NEO ACU-POINT TSUBO-KYU

A new type of indirect moxa. You are not required to push the moxa in the paper tube upward before use. You can quickly start using once you open up box. A gentle, and non-excessive heat penetrates the acu-point.

Type	pcs/box	Order Code
Gentle Heat	120	8009358
	600	8009359
Mild	120	8009360
	600	8009361



### NEO ACU-POINT MOXA CHOAN-DX

A new type of indirect moxa, just push the moxa in the paper tube upward, and apply directly to the skin.


Type	pcs/box	Order Code
Regular	120	8009305
	600	8009306
Mild	120	8009365
	600	8009366




# MOXA

## STICK-ON MOXA


**Soft Type**




**Light Type**



**Regular Type**



**Hard Type**



**CHOSEI-KYU MOXA**

Chosei-kyu is not burned directly on the skin surface but is raised above it so that its heat gives adequate stimulation. Heat regulating paper is enclosed for those with sensitive skin. Chosei-kyu's moxa burns evenly from top to bottom thus insuring a gradual heat build up that will not discomfort the skin. Each piece is self-adhesive and will hold firmly to skin. Simply remove from adhesive backing by pressing on reverse side with finger.

50 pcs. arranged on self-adhesive backing.

Type	pcs/box	Order Code
Soft*	100	8009370
	200	8009371
	1,000	8009372
Light	100	8009323
	200	8009324
	1,000	8009325
Regular	100	8009301
	200	8009302
	1,000	8009303
Hard	100	8009320
	200	8009321
	1,000	8009322

\*: without heat regulating paper

## HOW TO USE



- Remove individual moxa pieces by pushing on reverse side with your finger.
- Light moxa
- Apply to acupoints or painful areas. To regulate heat, apply regulating paper on skin first and apply moxa over paper.

### PRECAUTIONS

- The intensity of heat will vary depending on each person's physical condition and location of application. Please find your own most suitable level of warmth. If heat is too intense blisters may form, so remove from skin if too hot.
- Do not use moxa on the face.
- Do not use on wet skin as blisters may form. Dry skin thoroughly before application.
- The dried herbal compound will soften and its active ingredients melt at 40°C so please keep out of sunlight and store in a cool location. Keep out of reach of children.

Depending on constitution and skin sensitivity at point of application, the individual may be sensitive to the heat and blisters may form. In that case, please use the enclosed heat shield which will greatly decrease the heat intensity and make it difficult for blisters to form.

Caution: Do not touch flame directly to heat shield.

# MOXA

## DIRECT, LOOSE MOXA



**8009318:** 100 g/box  
**8009307:** 200 g/box



**8009354:** 100 g/box  
**8009351:** 200 g/box

### GOLD MOUNTAIN MOXA

Made from carefully selected superior quality mugwort. Super-pure, extra fine texture make this the highest quality moxa available. Excellent for direct moxibustion. Pure moxa is less likely to burn the skin because of its high concentration of essential oil. (Gold Mountain has a volume limit due to its limited gathering and production period.)

### FUKUJU MOXA

Medium grade moxa for everyday professional and home use. Made from good quality medium grade mugwort.

## DIRECT & NEEDLE MOXA



**8009312:** 100 g/box  
**8009313:** 200 g/box

### YOMEI MOXA

As a needle moxa it is of the highest quality with super cohesiveness for easy needle adhesion. As a direct moxa it is of standard quality, just below that of Fukuju and as a needle moxa is above Wakakusa.

## MINI MOXA SETS



### GOLD MOUNTAIN MOXA SET

Top quality moxa with 5 incense sticks for lighting. Contains 10 g of loose moxa. Excellent for home and clinic use and when travelling.

**8009315:** 10 g/box with 5 incense sticks



### FUKUJU MOXA SET

Standard quality moxa with 5 incense sticks for lighting. Contains 10 g of loose moxa

**8009326:** 10 g/box with 5 incense sticks

# MOXA

## NEEDLE MOXA



8009333



8009334

### HOT MOXA-BULK

For use in moxa needle heating. Bulk packing makes it economical. Two sizes are available,  $\phi$  12  $\times$  15 mm length or  $\phi$  10  $\times$  10 mm length each piece. Moxa pieces must be pierced with skewer or sharp object before use.

**8009333:** 500 pcs/box ( $\phi$  10  $\times$  10 mm length/box)

**8009334:** 500 pcs/box ( $\phi$  12  $\times$  15 mm length/box)



### WAKAKUSA MOXA

Wakakusa is the choice moxa for use in moxa needle heating (fire needle). Excellent adhesiveness with no stems. Its heat and easy lighting make it ideal for applying to acupuncture needles.

**8009309:** 300 g/box

**8009310:** 1 kg/box

**8009317:** 3 kg/box

## INDIRECT MOXA



### GOKUJO MOXA

Highest quality indirect moxa. Gokujo is ideal for burning on top of miso or sliced garlic and ginger in indirect moxibustion. Optimum quality and heating strength.

**8009319:** 300 g/bag

**8009367:** 5 kg/box



# MOXA

## INCENSE STICKS



### INCENSE STICKS FOR MOXA

Thick incense sticks ideal for use in Moxibustion. Handy and very easy to hold between fingers. Contains 10 bundles of incense sticks. 30 g/bundle including around 24 sticks.

**8009373:** Incense Sticks for Moxibustion, 10 bundles/box

## MOXA PRODUCTION PROCESS

**1**

Moxa is made from mugwort.

The underside of a mugwort leaf is covered with white hairs.

The collected white hairs become moxa.

**Picking**

Mugwort leaves are picked in summer.

**2**

Expose the leaves to sunlight and dry them in the sun for 1-2 days.

**Drying**

Apply hot air to the leaves for further drying. Reduce the water content of the leaves as much as possible.

**Grinding**

Grind the leaves finely in a millstone.

**3**

The leaves are then sifted through a traditional machine to separate the parts that will become moxa from the stems and other unwanted parts.

**Refining**

The moxa is mixed in a rotating bamboo basket for hours, and the powdery impurities are gradually removed.

Nagadoshi 長通し

Tomij 唐箕


**4**

The texture or purity of moxa is examined with the hand and eyes of experienced person in charge.

To maintain the consistency in the quality of final products, products from different production date are mixed.

# NON MAGNETIC PRESS PELLETS

## MAGRAIN ION PRESS PELLETS — Alloy Granules



### MAGRAIN ION PELLETS

Sakamura Magrain are used for both ear and body points. Magrain are available in 24K gold plated or silver color 1.2 mm alloy granules. Magrain is unique in that it provides an ionizing reaction which further enhances the compression stimulation to the point. Each pellet comes with a new high quality skin (tan) colored or clear (transparent) round  $\phi$  7 (0.28" dia.) hypo-allergenic band-aid tape.

300 pellets/pack

Type	Tape	Pack	Order Code
Silver	with clear	with pack	8010106
	with clear	without pack	8010105
	with tan color	with pack	8010102
	with tan color	without pack	8010101
Gold	with clear	with pack	8010108
	with clear	without pack	8010107
	with tan color	with pack	8010104
	with tan color	without pack	8010103



### MAGRAIN-5 MULTI-PELLETS

Each 12 × 12 mm (1/2" × 1/2") tan colored band-aid tape plaster features a geometrically space cluster of five  $\phi$  1.2 alloy granules. The clusters act to stimulate a wider treatment area.

24K gold plated Magrain-5 also available.

50 plasters/pack

Type	Pack	Order Code
Silver	with pack	8010110
	without pack	8010109
Gold	with pack	8010123
	without pack	8010124

## HOW TO USE PRESS PELLETS



Magrain and Accu-Patch (see next page) are commonly used for both ear and body points. The small  $\phi$  1.2 stainless steel, silver or gold pellets produce a long acting, continuous compression effect on the point. Gold is said to tonify or energize; Silver sedates or calms; and Stainless Steel balances or harmonizes.

Generally, pellets should be changed every 2 to 3 days to achieve the maximum effect (usually the body will adapt to the pressure after 2 to 3 days). In regards to conical pellets, it is recommended to be changed every 6 hours after use. Before applying, the skin should be thoroughly cleaned with alcohol and the tape applied with a good seal.

The patient may be reminded to gently press or stimulate these pellets occasionally during the day or as otherwise directed.



### MAGRAIN CONICAL PELLETS


The conical, 2 mm height stainless steel pellets provide more stimulation to the points, and you can expect immediate effect. Each pellet comes with high quality tan colored round shape  $\phi$  12 (0.47" dia.) plaster. It is applied for both ear and body points likewise the standard pellets from Magrain Gold and Silver types. Due to its strong stimulation, it is recommended to change pellets every 6 hours.

50 pellets/pack

Type	Pack	Order Code
Conical	with pack	8010152
	without pack	8010151

# NON MAGNETIC PRESS PELLETS

## ACCU-PATCH PRESS PELLETS - Stainless Steel Granules




### ACCU-PATCH STAINLESS STEEL PRESS PELLETS

New improved stainless steel press pellets for the ear and body. Similar to Magrain except the  $\phi$  1.2 stimulating press ball is made from high quality surgical stainless steel. Each pellet comes with a new high quality tan (skin) colored or clear (transparent) round  $\phi$  7 (0.28" dia.) hypo-allergenic band-aid tape.

300 pellets/pack

**8010119:** Accu-Patch with clear tape

**8010120:** Accu-Patch with tan color tape




### ACCU-PATCH GOLD PLATED PRESS PELLETS

24K gold plated press pellets.

300 pellets/pack

**8010121:** Accu-Patch Gold Plated with clear tape

**8010122:** Accu-Patch Gold Plated with tan color tape




### ACCU-PATCH CONICAL PELLETS

The conical shape of the pellets provide more stimulation to the points, and you can expect immediate effect. The pellet is made of stainless steel, 2 mm high, and the plaster is  $\phi$  12 (0.47" dia.).

50 pellets/pack

**8010162:** Accu-Patch Conical Pellets with tan color tape

## ACCU-PATCH TITANIUM PELLETS-High Purity Titanium Granules



### ACCU-PATCH TITANIUM PRESS PELLETS


In the medical field, titanium's high level of bio-compatibility has led to its use in many instruments, such as artificial bones, joints, etc. However, the production of smooth, high-purity pellets of the appropriate size has long been a challenge to the industry. After years of research and testing, a new method has been developed in which titanium pellets are produced by melting titanium with high temperature plasma in a clean, inert gas. This special method enables high purity titanium to be produced for use in the acupuncture pellets. These pellets can be used in treatment without any fear of disintegration caused by impurities or oxidization. The pellets are hypoallergenic and are free from unwanted side effects. This Titanium Accu-Patch pellets have been approved by the Japanese Ministry of Health and Welfare in recognition of their outstanding safety and stability.

Tape	pellets/pack	Order Code
clear	200	8010129
	300	8010127
tan color	200	8010130
	300	8010128



# MAGNETIC PRESS PELLETS

## ACCU-BAND 800 - FERRITE MAGNETS




**ACCU-BAND** emits a constant pressure and a concentrated magnetic flux of 80 mT or 800 gauss, thus giving continuous effect to shoulders, arms, legs or wherever ACCU-BAND is applied. Adheres even in water and will not rust. The ferrite magnet features a mixed north-south pole structure without clearly defined pole surfaces.


**012486:** Non-plated, 2 sheets (24 pcs)/box

**012488:** Gold-plated, 2 sheets (24 pcs)/box

**012485:** Non-plated, bulk, 12 pcs/sheet





## ACCU-BAND HIGH STRENGTH REC MAGNETS



### ACCU-BAND 6000

New 600 mT or 6,000 gauss (Br) high strength magnets. 24K gold plated  $\phi$  4.0 x 1.6 mm thick permanent magnets with  $\phi$  19 (3/4" dia.) band-aids. The smooth surface of the north pole faces upward, and the underside, representing the south pole, is indicated by a dent.

1 sheet (12 pcs)/box




### ACCU-BAND 9000

New 900 mT or 9,000 gauss (Br) high strength magnets. 24K gold plated  $\phi$  5.0 x 1.6 mm thick permanent magnets with  $\phi$  19 (3/4" dia.) band-aids. A little larger in diameter than ACCU-BAND 6000. The smooth surface of the north pole faces upward, and the underside, representing the south pole, is indicated by a dent.

1 sheet (12 pcs)/box

## REPLACEMENT PLASTERS



A hole in the plaster backing paper allows the press needle to be easily and accurately placed on the plaster.

012494 are intended as replacements for the MAG-PRESS needle or any kind of small magnet.

### 012494:

Round shape,  $\phi$  19

12 plasters/sheet

10 sheets (120 plasters)/pack

### 271064:

Round shape,  $\phi$  19

12 plasters/sheet

100 sheets (1,200 plasters)/pack

# ACUPUNCTURE DEVICES

## ES-160 6-CHANNEL UNIT with RYODORAKU FUNCTION



ES-160 complete with standard accessories, including

<b>011380</b>	Search/Stimulation Probe	<b>011328</b>	Electrode Cable (Blue)
<b>011324</b>	Electrode Cable (Black)	<b>011329</b>	Electrode Cable (Grey)
<b>011325</b>	Electrode Cable (Brown)	<b>B060999</b>	Exclusive Clips (Black), 6x
<b>011326</b>	Electrode Cable (Red)	<b>B061000</b>	Exclusive Clips (Red), 6x
<b>011327</b>	Electrode Cable (Green)		

NOTE: Batteries are not included.

### for Ryodoraku treatment, Electroacupuncture, and Acupuncture point search/stimulation

This is an ideal model, designed to respond to the needs of practitioners in the field, including such features as Ryodoraku measurements, preset programs, and free programmable memories, to name but a few.

- Allows measurement of Ryodo Points and storage of the measurement data.
- Allows simultaneous stimulation with twelve (12) acupuncture needles (6 channels).
- Provides eight (8) different stimulation modes, including three (3) preset program modes.
- Saves parameters of the previous treatment and incorporates sixteen (16) free programmable memory settings.
- Enables user to find and stimulate acupuncture points.

### About skin conductivity measurement:

ES-160 is unique in the scientific application of Oriental Medicine. Oriental Medicine recognizes that an electro-potential relationship exists between the acu-point on the skin and the internal organs.

### Specifications

Power source: Four(4) R14 battery  
Ryodoraku measurement range: 0–200  
Voltage setting: 6, 12 or 21 volt  
Number of channels: 6, independently controllable  
Pulse shape: Symmetric, bi-phasic square pulse  
Frequency: 0.3–500 Hz  
Phase duration: 50–400  $\mu$ s, adjustable  
Amplitude:

- Stimulation by 'search and stimulate' probe: 0–48 mA (peak)
- Electro-acupuncture:
  - High: 0–32 mA, Low: 0–16 mA (peak)

Times: 0–60 minutes

Output modes:

- (1) Constant (2) Burst (3) Surge (4) Fast + Slow  
(5) Sweep (6) Random prog. 1 (7) Random prog. 2  
(8) Random prog. 3

Number of memory setting:

- Memory for Ryodoraku measurement: 24
- Memory settings for previous treatment parameters: 5, one for each mode (random modes excluded)
- Free programmable memory settings: 5

Dimensions: 239(W) x 174(L) x 41(H) mm

Weight: Approx. 600 g (without batteries)

# ACUPUNCTURE DEVICES

## ES-130 3-CHANNEL PALM-SIZED ELECTROACUPUNCTURE



### High Performance and Versatility

- Low/High intensity settings
- Frequency selector offers a complete frequency range (1–500 Hz)
- Suitable for needles

### Unique Wave Form

- Unique wave pulse minimizes excessive charge build up in the tissue and damaging the nerve

### Safe, Reliable, and Easy to Use

- Each of the 3 channels are independently controlled
- Special dial controls allow easy use yet safe operation
- Built in battery life check with indicator light



260182



ES-130 complete with standard accessories, including

**260182** Carrying Bag

**B180432** Electrode Cable with Clips (white)

**B180433** Electrode Cable with Clips (yellow)

**B180434** Electrode Cable with Clips (green)

NOTE: A battery isn't included.



013031



013032



013033

013034

Optional accessories

**013031** Belt Clip with Screws

**013032** Electrode Cable

**013033** Exclusive Clips (Black)

**013034** Exclusive Clips (Red)

### Specifications

Power supply: 1 × DC 9 V battery

Frequency: 1–500 Hz, variable

Maximum output voltage: 30 V (at 500  $\Omega$  load)

Current amplitude (max.): 14 mArms (at 500  $\Omega$  load)

Phase duration: 100  $\mu$ s

Number of channel: 3

Others: Battery checker Pulse rate indication

Safety class according to IEC 60601-1: Internally powered equipment, Type BF

Dimensions: 63(W) × 27(D) × 96(H) mm

Weight: 120 g (without a battery)

# KINESIO ELASTIC TAPE

## OLYMPIA KINESIOLOGY TAPE



### What is Kinesiology Tape?

Kinesiology Tape is used to aid muscle or joint healing as well as give support and stability without restrictions to movement or circulation. Best for injured muscles or joints, swelling and pain relief, all while staying active!

It is water-resistant, breathable and as flexible as human skin, making it capable of staying on up to several days. When applied correctly, the taped area forms creases expanding the space between skin and muscle to release pressure and improve fluid circulation.



- Pre-cut and Roll types are available.
- Eight color choices (Beige, Blue, Pink, Black, Yellow, Green, Orange and Red).
- Various tape width (2.5 cm, 3.75 cm, 5 cm and 7.5 cm).
- Compact-sized carry box.
- Optimized elasticity gives stronger and gentle support.
- Improved fabric, extra comfort, gentle for sensitive skin.
- New synthetic engineered Microfiber technology. Reflect under the lights.
- Water-proof fabric, instant dry.
- Medical grade, heat sensitive acrylic adhesive. (Meet ISO10993 requirement)



Length	Type	Color								roll(s)/box	Min. Order
		Beige	Blue	Pink	Black	Yellow	Green	Orange	Red		
5 cm x 5 m	Roll	8011401	8011402	8011403	8011404	8011405	8011406	8011407	8011408	1	10 boxes
5 cm x 25 cm/pc	Pre-cut	8011411	8011412	8011413	8011414	8011415	8011416	8011417	8011418	1	10 boxes
7.5 cm x 5 m	Roll	8011421	—	—	—	—	—	—	—	4	10 boxes
5 cm x 5 m	Roll	8011422	—	—	—	—	—	—	—	6	10 boxes
3.75 cm x 5 m	Roll	8011423	—	—	—	—	—	—	—	8	10 boxes
2.5 cm x 5 m	Roll	8011424	—	—	—	—	—	—	—	12	10 boxes
5 cm x 31 m	Roll	8011428	8011432	8011433	8011434	8011435	8011436	8011437	8011438	1	10 boxes



# KINESIO ELASTIC TAPE

## SARASA KINESIOLOGY TAPE



SARASA tape treats a wide variety of conditions due to its ability to alleviate pain, reduce inflammation, relax the muscles, enhance performance and facilitate while supporting the muscle in movement.

- 100% cotton tape material, latex free
- 30–40% elasticity from resting length
- Medical grade, thermoplastic Skin-Care adhesive
- Water repellent
- Allows skin to breathe

Length	Color				roll(s)/box	Min. Order
	Beige	Black	Blue	Red		
5 cm x 5 m	8011310	8011308	8011307	8011306	1	100 boxes
7.5 cm x 5 m	8011316	—	—	—	4	20 boxes
5 cm x 5 m	8011311	8011305	8011304	8011303	6	20 boxes
3.75 cm x 5 m	8011309	—	—	—	8	20 boxes
2.5 cm x 5 m	8011314	—	—	—	12	20 boxes
5 cm x 30 m	8011315	—	—	—	1	20 boxes

## SARASA COTTON SUPPORT TAPE



COTTON SUPPORT TAPE can be used to fix joints and muscles.

- Special made highly breathable fabric.
- Good ventilation to prevent eruption.
- Easy to storage, easy to cut.
- Long lasting adhesiveness.




Length	roll(s)/box	Order Qty.	Order Code
3.8 cm x 13.7 m	8	20 boxes or more	8011321
		1–19 box	8011322
5 cm x 13.7 m	6	20 boxes or more	8011323
		1–19 box	8011324

## SARASA UNDERWRAP



UNDERWRAP can be used to protect skins. It should be applied on the skin before applying COTTON SUPPORT TAPE or KINESIOLOGY TAPE.

- Good ventilation to prevent eruption.
- Enough flexibility, elasticity and durability.
- Firmly protect skin even it's quite thin.



Length	roll(s)/box	Order Qty.	Order Code
7 cm x 27.4 m	12	16 boxes or more	8011325
		1–15 box	8011326

# KINESIO ELASTIC TAPE

## SARASA NEVANON SCISSORS



8011317



8011318




8011319

NEVANON SCISSORS is specialized to cut COTTON SUPPORT TAPE, UNDERWRAP and KINESIOLOGY TAPE easily.

- Fluorine resin coating on the inside of blade to cut tapes easily.
- 3 size choices (small, medium and large).

**8011317:** Small  
**8011318:** Medium  
**8011319:** Large

- **Long-lasting smooth cutting is available with the fluorine resin coating.**



GENERAL CATALOGUE  
for **ACUPUNCTURE**

**ITO**co.,LTD.

3-1-8 Sakae-cho, Kawaguchi-shi, Saitama 332-0017, Japan  
TEL: 81-48-254-1031 FAX: 81-48-254-1033  
URL: <https://www.itocoltd.com/> E-Mail: [itocoltd@itolator.co.jp](mailto:itocoltd@itolator.co.jp)

\*Designs and specifications are subject to change without notice.

No.162 C0c000256-2505 Z2505DATA